



novigenix

LITOSeek™

Immunotherapy Response Prediction

Precise & Predictive Patient Profiling

Introduction

Traditional tumor biopsy techniques, while informative, face limitations due to tumor heterogeneity and accessibility challenges at various metastatic sites throughout a patient's treatment journey. Advances in next-generation sequencing technology (NGS) have enabled the development of non-invasive methods for blood-based biomarker detection, commonly referred to as a liquid biopsy. Novigenix's pioneering liquid biopsy technology platform LITOSseek™ (Liquid Immuno-TranscriptOmics), leverages whole-blood RNA-sequencing to capture patient-specific transcriptomic profiles. By integrating next-generation sequencing and AI technologies, LITOSseek™ excels in discovering robust blood-based RNA biomarkers and creating precise, specialized predictive models. This application note highlights the utility and benefits of LITOSseek™ in clinical trials including prediction of treatment response to immunotherapy with durable therapeutic benefit.





Figure 1: Schematic overview of the LITOSseek™ platform, illustrating workflow from blood sample collection to the construction of precise predictive models.


Prediction of durable therapeutic benefit

Bladder cancer patients are frequently treated with anti-PD-1/PD-L1 immunotherapy as a second line treatment option. However, only a small percentage (15-24%) experience a durable benefit. LITOSeek™ was employed to analyze whole blood samples from 89 metastatic bladder cancer patients, both before and during immunotherapy. The patients were divided into three independent cohorts for discovery, confirmation, and validation, with the former aiding in biomarker identification and model training, and the latter for independent validation of predictive models.


Blood based longitudinal immunotherapy response prediction



A total of 89 metastatic bladder cancer patients underwent immunotherapy over four cycles, and clinical responses were evaluated six months after therapy initiation. Whole-blood PAXgene® Blood RNA Tubes were collected before treatment and on-treatment (either after the first or second cycle), and RNAseq data was generated for analysis to identify gene expression signatures predictive of durable therapeutic benefits to immunotherapy and characterization of immune modulation in patients with metastatic bladder cancer.




The LITOSseek™ NoviScore has been optimized to accurately pinpoint the most significant Differentially Expressed Genes (DEGs) as shown in this representative volcano plot of NoviScore for discovery of key biomarkers indicative of immune changes associated with effective response to therapy.




The DEGs identified with LITOSseek™ provide deep insights into the immune system's dynamic response to immunotherapy. Moreover, the results are consistent with multiple other studies concluding that response to immunotherapy is driven by systemic cell responses, with a major role played by T cells which can be predictive of clinical benefit (Fairfax et al., Nature Medicine 2020; Valpione et al, Nature Cancer 2020; Hiam-Galvez, K.J. et al., Nat Rev Cancer, 2021).

NoviContext uses multiple knowledge-based bioinformatic approaches, such as protein-interaction network analysis and biological pathway analysis to characterize central pathways related to early response to immunotherapy. Contextualization in this cohort has identified key pathways of antigen processing and presentation, T-cell mediated immunity, and cell cycle in patients with durable therapeutic benefit.





LITOSseek™ AI-enabled modeling precisely predicts response to ICI

LITOSseek™ NoviSig feature reduction machine-learning algorithms identified a set of 10 genes related to relevant group comparisons and central biological pathways for downstream modeling. Predictive models developed based on these 10 early-response biomarkers were trained with NoviBoost algorithms, and further validated on an independent cohort. The Kaplan-Meier curve illustrates the predictive power of LITOSseek™ in stratification of patients into responders and non-responders to immunotherapy:



These results highlight the platform's capabilities in identifying highly accurate RNA biomarkers of therapeutic response where other traditional biomarkers such as PD-L1 Score and Tumor Mutation Burden fail to have predictive power:



Conclusions

The LITOSseek™ platform has demonstrated a broad range of applicability in precision oncology, offering Precise & Predictive Patient Profiling such as prediction of therapeutic response to immunotherapy, thereby positioning it as a transformative tool in the landscape of novel therapy development. Other platform capabilities available for clinical development of novel therapies include:



**Immuno
pharmaco-dynamics
(ImmunoPD)
response markers**



**Predictive
biomarkers of
response and
resistance to therapy**



**Mechanisms
of action and
mechanisms of
resistance to therapy**



**Discovery of safety
biomarkers**



**Digital cytometry and cellular
deconvolution overcoming
logistical hurdles of flow
cytometry and CyTOF**



Contact us

Route de la Corniche 3
Phenyl Building
CH-1066 Epalinges
Switzerland

✉ info@novigenix.com

Follow Us

 [novigenix](#)

novigenix.com